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
Peer-reviewed article: 2021-09-01

Perioral Aerosol Sequestration Suction Device Effectively Reduces Biological Cross-Contamination in Dental Procedures

Victor Usón¹, María Laura González², Vicente Lozano-de Luaces³, María Cristina Martínez⁴

Abstract

Objective: The infection risk during dental procedures is a common concern for dental professionals which has increased due to coronavirus severe acute respiratory syndrome coronavirus 2 pandemic. The development of devices to specifically mitigate cross-contamination by droplet/aerosol is crucial to stop infection transmission. The objective of this study is to assess the effectiveness of a perioral suction device (OBF) to reduce biological contamination spread during dental procedures.

Materials and Methods: Forty patients were randomized 1:1 to a standard professional dental hygiene treatment with OBF and without. Aerosolized hydrophobic (HPL) fluorescence assay was used to evaluate the spread of potential contaminants. The total number of visible light units (VLU) from key dental operator locations: operator's face-mouth, back of the surgical operator's gloves, patient's safety goggles, and environmental table were measured. Percentage contamination reduction between control and OBF were compared.

Statistical Analysis: Primary outcome, total VLU, was analyzed by comparing the means of logged data, using a two-sided two-sample t-test. Secondary outcomes are VLU of logged data for the different locations were analyzed in the same way. Proportion of patients from whom different locations reported events (clean, acceptable, and failed) were analyzed by using Fisher's exact test.

Results: For the whole dental environment, VLU reduction (>150 units) achieved with OBF was 98.4% (97.4%-99.5%). By dental operator location the reduction in VLU

Acceso al estudio. Fuente NCBI: National Center for Biotechnology Information